# Catalogue PolyDiagnost

Version: 2.7

With publication of this catalogue older versions loose their validation

POLYDIAGNOST Am Söldnermoos 17 85399 Hallbergmoos Germany

Phone: +49 811 99873380 Fax: +49 811 99873389 info@polydiagnost.com www.polydiagnost.com

09.04.2015



## Table of contents

Optics	2
Accessories for Optics	4
Combination table semirigid Endoscopes	7
Semirigid Endoscopy	8
Combination table PolyScope	18
Flexible Endoscopy	19
Electronics	
Acessories for Electronics	24
Spare Parts	26



Optics			
Order No.	Article Description	Additional Description	Picture

PD-DS-1080	semirigid optic modular 30,000 Pixel / 70°	the interdisciplinary semirigid endoscope one optic for every application highest resolution examination on the highest level no decontamination issues – no lumen! 30,000 pixel / 70° FOV outer diameter: 1.3 mm standard light connection: ACMI/Storz/Wolf overall length: 181mm only on request	
PD-DS-1083	semirigid optic modular 10,000 Pixel / 120°	the interdisciplinary semirigid endoscope one optic for every application high resolution at low diameter, 120° FOV  no decontamination issues – no lumen! 10,000 pixel / 120° FOV outer diameter: 0.9 mm standard light connection: ACMI/Storz/Wolf overall length: 181mm	<u>O</u>
PD-DS-1084	semirigid optic modular 6,000 Pixel / 70°	the interdisciplinary semirigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues – no lumen! 6,000 pixel / 70° FOV outer diameter: 0.55 mm standard light connection: ACMI/Storz/Wolf overall length: 181mm	Q
PD-DS-1086	semirigid optic modular 3,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application ultrathin, with finest resolution especially for diagnostic in thinnest channels no decontamination issues – no lumen ! 3,000 pixel / 70° FOV outer diameter: 0.45 mm standard light connection: ACMI/Storz/Wolf overall length: 181mm	Q."
PD-PN-0080	semirigid optic modular PCNL 10,000 pixel / 120°	semirigid endoscope special design for PCNL high resolution at low diameter, 120° FOV no decontamination issues – no lumen ! 10,000 pixel / 120° FOV outer diameter: 0.9 mm standard light connection: ACMI/Storz/Wolf overall length: 310mm	Q_



Optics			
Order No.	Article Description	Additional Description	Picture

PD-PS-0095	flexible optic PolyScope 120° view / 300cm	PolyScope- the flexible optic system light fibres are integrated in the disposable catheter full light in every adoption no repairs needed no decontamination issues due to disposable catheters 10,000 Pixel, 120° FOV, 300cm length suitable catheter length: 185cm	0
PD-PS-0096	flexible optic PolyScope 120° view / 180cm	PolyScope- the flexible optic system light fibres are integrated in the disposable catheter full light in every adoption no repairs needed no decontamination issues due to disposable catheters 10,000 Pixel, 120° FOV, 180cm length suitable catheter length: 70cm	
PD-PS-0097	flexible optic PolyScope 120° view / 150cm	PolyScope- the flexible optic system light fibres are integrated in the disposable catheter full light in every adoption no repairs needed no decontamination issues due to disposable catheters 10,000 Pixel, 120° FOV, 150cm length suitable catheter length: 30 and 42cm	Q



PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm optimized for use with: 10,000 and 30,000 pixel optic	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm optimized for use with: 6,000 and 3,000 pixel optic	=mmp
PD-AC-0092	3-joint-arm friction	with knob for use with semirigid endoscopes  - fixes the ocular which is connected to the optic body and the light cable  - can be fixed in every position and serves as a"third" hand keeps all accessories in the unsterile area	
PD-AC-0093	3-joint-arm friction PolyScope	with knob for use with flexible endoscopes  - fixes the ocular which is connected to the optic body and the adaptor seperate light  - can be fixed in every position and serves as a "third" hand  - keeps all accessories in the unsterile area	
PD-AC-0091	3-joint-arm friction short	similar to PD-AC-0092 with short camera support recommended for sales agents	
PD-AC-0095	3-joint-arm friction PolyScope short	similar to PD-AC-0093 with short camera support recommended for sales agents	301
PD-AC-1012	SteriTray large	SteriTray for PCNL optic with individual holders for the endoscope size: 250x480x60mm	



### Catalogue

# Accessories for Optics Order No. Article Description Additional Description Picture

PD-AC-1112	cover plate large	cover plate for SteriTray PD-AC-1012 recommended for the secure transport	
PD-AC-1013	SteriTray small	SteriTray for semirigid optics (except PD-PN-0080) with individual holders for the endoscope size: 240x250x60mm	
PD-AC-1113	cover plate small	cover plate for SteriTray PD-AC-1013 recommended for the secure transport	





## The modular semirigid system

- no decontamination issues
- cost- effective
- interdisciplinary




For possible combinations and fields of applications see combination table next page



### **Combination table semirigid Endoscopes**

The interdisciplinary semirigid endoscopy system- one optic for every application!

Without any decontamination issues due to the sterile, disposable PolyShaft.



Please notice: errors excepted; subject and prices to change without notice; similar illustrations



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	

PD-AC-0050	Solex-Nipple- Expander 1.15mm	soft-lumen-expander suitable for PolyShaft 0.90 mm essential for milk duct endoscopy; alleviates the insertion of the endoscope in the milk duct please order separately: obturator PD-AC-0055	0 10
PD-AC-0051	Solex-Nipple- Expander 0.95mm	soft-lumen-expander suitable for PolyShaft 0.65 and 0.75 mm essential for milk duct endoscopy; alleviates the insertion of the endoscope in the milk duct please order separately: obturator PD-AC-0055	
PD-AC-0052	Solex-Nipple- Expander 1.3mm	soft-lumen-expander suitable for for PolyShaft 1.15 mm and viewing forceps essential for milk duct endoscopy; alleviates the insertion of the endoscope in the milk duct please order separately: obturator PD-AC-0055	3 = = =
PD-AC-0055	obturator for Solex	required for all Solex-nipple-expanders disposable / sterile packed packing unit: 50pcs	
PD-CA-0113	Tuohy Borst adaptor	disposable / sterile packed	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0203	Y-adaptor	disposable / sterile packed	



Semirigid Endoscopy			
Order No.	Article Description	Additional Description	Picture

PD-CA-0606	check valve (Schanuppsi)	required for drill PD-DS-1620,PD-DS-1630 and PD-DS-1640 and also recommended for baskets and viewing forceps  disposable / sterile packed	
PD-CA-0609	Luer Lock cap set	set 5 pieces disposable / sterile packed	-
PD-CA-0610	Acha 6000 Pixel	for use of semirigid optic without PolyShaft suitable for 6.000 pixel optic diameter: 0.59mm LuerLock	
PD-CA-0611	Acha 10000 Pixel	for use of semirigid optic without PolyShaft for 10.000 pixel optic diameter: 0.95mm LuerLock	
PD-CA-0612	Acha 30000 Pixel	for use of semirigid optic without PolyShaft for 30.000 pixel optic diameter: 1.33mm LuerLock only on request	
PD-DS-0500	lavage pipe 0.5mm	Lavagepipe incl. check valve outer diameter: 0.5mm targeted lavage with the semirigid modular system- two channels in one PolyShaft disposable / sterile packed	
		only on request	
PD-DS-0505	irrigation-/suctionpipe 0.8mm	irrigation-/suction pipe incl. check valve outer diameter: 0.8mm targeted irrigation and suction with the semirigid modular system- two channels in one PolyShaft	-
		disposable / sterile packed	
		only on request	



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	

PD-DS-1000	PolyShaft straight 0.65 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 0.55/0.65mm working length: 90mm autoclavable packing unit: 10pcs	
PD-DS-1005	PolyShaft straight 0.75 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 0.65/0.75mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-DS-1010	PolyShaft straight 0.90 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 0.70/0.90mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-DS-1015	PolyShaft straight 1.15 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 1.05/1.15mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-DS-1020	PolyShaft straight 1.60 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 1.40/1.60mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-DS-1025	PolyShaft straight 2.00 / 90 mm	special design for easy insertion of the optic inner/outer diameter: 1.80/2.00mm Working length: 90mm autoclavable packing unit: 10pcs	
PD-DS-1030	PolyShaft straight 0.75 / 50 mm	special design for easy insertion of the optic inner/outer diameter: 0.65/0.75mm Working length: 50mm disposable / sterile packed packing unit: 100pcs	



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	

PD-DS-1035	PolyShaft straight 0.90 / 50 mm	special design for easy insertion of the optic inner/outer diameter: 0.70/0.90mm Working length: 50mm disposable / sterile packed packing unit: 100pcs	00
PD-DS-1038	PolyShaft straight 1.10/ 50 mm	special design for easy insertion of the optic inner/outer diameter: 1.00/1.10mm Working length: 50mm disposable / sterile packed packing unit: 100pcs	-
PD-DS-1185	trocar sharp	for PolyShaft 16G (PD-DS-1042) set of 10 pieces autoclavable	
PD-DS-1195	obturator blunt	for PolyShaft 16G (PD-DS-1042) set of 10 pieces autoclavable	
PD-DS-1210	optic shifter	length adjustment for all semirigid optics and handles length 26mm autoclavable	
PD-DS-1220	optic shifter ophthalmology, short	length adjustment for semirigid optics and short handle length 48mm autoclavable	
PD-DS-1315	handle middle	suitable for PolyShaft 50mm 3 female Luer Lock connections for irrigation and instrumentation  autoclavable	
PD-DS-1320	handle short	suitable for PolyShaft 90 mm 3 female Luer Lock connections for irrigation and instrumentation  autoclavable	

www.polydiagnost.com



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	

PD-DS-1620	spiral drill 50	spiral drill for handle middle with PolyShaft 50 mm diameter: 0.38 mm / length over all: 192 mm only to use with check valve PD-CA-0606 (1 incl.)	
PD-DS-1630	spiral drill 90	spiral drill for handle short with PolyShaft 90	Mallo Medical Control of the Control
NEW	·	mm diameter: 0.38 mm / length over all: 199 mm only to use with check valve PD-CA-0606 (1 incl.)	
		autoclavable	
PD-DS-1640	spiral drill 50 length 158 mm	spiral drill for handle short and Shifter 48 mm with PolyShaft 50 mm  - ophthalmology  diameter: 0.38 mm /length over all: 158 mm only to use with PD-CA-0606 (1 included)	
		autoclavable	
PD-OI-0001	dilator (according to Wilder)	dilator according to Wilder  autoclavable	
PD-OI-0010	Bowman probe 00/0	buttoned shape Fig. 00/0 diameter: 0.7/0.8 mm  autoclavable	
PD-OI-0011	Bowman probe 0/1	buttoned shape Fig. 0/1 diameter: 0.8/0.9 mm  autoclavable	
PD-OI-0012	Bowman probe 1/2	cylindrical shape Fig. 1/2 diameter: 0.9/1.1 mm  autoclavable	



PD-OI-0013	Bowman probe 2/3	cylindrical shape Fig. 2/3 diameter: 1.1/1.3 mm  autoclavable	
PD-PN-0020	10.5 F working shaft	working shaft 10.5F for the PCNL system suitable for the lithoclast® master probe inner/outer diameter 8.5F/10.5F including Dilator  autoclavable	
PD-PN-0030	2-port-adaptor for PCNL	2-port-adaptor for PCNL working shaft 10.5F with irrigation port  autoclavable	
PD-PN-1008	MicroPerc® Set complete	content: 3-port with female Luer-Lock, Luer-Lock cap, Tuohy Borst adaptor, perfusion tube, 3-way stopcock and 8F working shaft, Allseeing-needle, working shaft 4.8 F working shaft 8F for the PCNL system suitable for the lithoclast® probe inner/outer diameter 7F/8F including Dilator  set of 5 pieces	ELIOAGHOIT
PD-PS-0124	lasershifter 1000	laser- and instrument shifter for outer diameter 500-1000µ  only on request	
PD-PS-0125	lasershifter 500	laser- and instrument shifter for outer diameter 200-500µ	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush for optic shifter PD-DS-1210 set of 5 pieces	0
PD-TD-0030	cleaning brush handle set of 5	cleaning brush for all handles set of 5 pieces	0

Please notice: errors excepted; subject and prices to change without notice; similar illustrations



PD-TD-0040	cleaning brushes for 2-port-adaptor for PCNL	cleaning brush for PCNL 2-port-adaptor set of 2x5 brushes	0
PD-TD-0060	cleaning brush for 10,5 F workingshaft	cleaning brush for PCNL working shaft 10.5F set of 2x5 brushes	
PD-TI-2106	handle for modular basket with chuck	handle with chuck suitable for the modular baskets PD-TI-3503 and PD-TI-3504  autoclavable	
PD-TI-3503	3 wire basket	3 wire basket without handle for semirigid endoscopy diameter: 0.38mm / length: 350mm disposable / sterile packed please order separately: PD-TI-2106	
PD-TI-3504	4 wire basket	4 wire basket without handle for semirigid endoscopy diameter: 0.38mm / length: 350mm disposable / sterile packed please order separately: PD-TI-2106	
PD-TI-5110	grasping forceps semirigid	grasping forceps for semirigid endoscopy with scissor handle and toothed spoon diameter: 0.8 mm, length: 35 cm  autoclavable	<b>√</b>
PD-TI-8101	biopsy forceps semirigid	biopsy forceps for semirigid endoscopy with scissor handle and oval spoon diameter: 0.8 mm, length: 35 cm  autoclavable	90



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	

PD-TI-9901	oil spray	oil spray for maintenance suitable for all forceps	
PD-ZS-0913	Dilatation Balloon 30mm	dilatation ballon for semirigid endoscopy balloon length: 30 mm  disposable / sterile packed	
PD-ZS-9142	Saliva Drain Type I 4.5F 20	salive drain type I for postoperative splinting of the saliva channel diameter: 4.5 F, length: 20 mm disposable / sterile packed	
PD-ZS-9145	Saliva Drain Type I 4.5F	salive drain type I for postoperative splinting of the saliva channel diameter: 4.5 F, length: 40 mm disposable / sterile packed	
PD-ZS-9160	Saliva Drain Type I 6F	salive drain type I for postoperative splinting of the saliva channel diameter: 6 F, length: 40 mm disposable / sterile packed	
PD-ZS-9162	Saliva Drain Type I 6F 20	salive drain type I for postoperative splinting of the saliva channel diameter: 6 F, length: 20 mm disposable / sterile packed	
PD-ZS-9190	Saliva Drain Type I 9F	salive drain type I for postoperative splinting of the saliva channel diameter: 9 F, length: 40 mm disposable / sterile packed	



Semirigid Endoscopy				
Order No.	Article Description	Additional Description	Picture	


PD-ZS-9192	Saliva Drain Type I 9F 20	salive drain type I for postoperative splinting of the saliva channel diameter: 9 F, length: 20 mm  disposable / sterile packed	
PD-ZS-9245	Saliva Drain Type II 4.5F	salive drain type II for postoperative splinting of the saliva channel diameter: 4.5 F, length: 40 mm disposable / sterile packed	
PD-ZS-9260	Saliva Drain Type II 6F	salive drain type II for postoperative splinting of the saliva channel diameter: 6 F, length: 40 mm disposable / sterile packed	
PD-ZS-9290	Saliva Drain Type II 9F	salive drain type II for postoperative splinting of the saliva channel diameter: 9 F, length: 40 mm disposable / sterile packed	
PD-ZS-9345	Saliva Drain Type III 4.5F	salive drain type III for postoperative splinting of the saliva channel diameter: 4.5 F, length: 120 mm can be shortened from the distal end disposable / sterile packed	
PD-ZS-9360	Saliva Drain Type III 6F	salive drain type III for postoperative splinting of the saliva channel diameter: 6 F, length: 120 mm can be shortened from the distal end disposable / sterile packed	
PD-ZS-9390	Saliva Drain Type III 9F	salive drain type III for postoperative splinting of the saliva channel diameter: 9 F, length: 120 mm can be shortened from the distal end disposable / sterile packed	





# PolyScope – the flexible modular system


- no decontamination issues
- no more repairs due to broken light fibres
- interdisciplinary



For possible combinations and fields of applications see combination table next page



### **Combination table PolyScope**





# Flexible Endoscopy Order No. Article Description Additional Description Picture

PD-CA-0108	Y-adapter with haemostasis valve	disposable / sterile packed now without obturator	W/W
PD-CA-0113	Tuohy Borst adaptor	disposable / sterile packed	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0606	check valve (Schanuppsi)	required for drill PD-DS-1620, PD-DS-1630 and PD-DS-1640 and also recommended for baskets and viewing forceps  disposable / sterile packed	
PD-LC-9551	adaptor separate light	connects the PolyScope catheter to the light source with pin for easy mounting at the 3- jointarm standard light connection: ACMI	
PD-PS-0120	optic combishifter	length adjustment for the PolyScope	
PD-PS-0124	lasershifter 1000	laser- and instrument shifter for outer diameter 500-1000µ  only on request	
PD-PS-0125	lasershifter 500	laser- and instrument shifter for outer diameter 200-500µ	



# Flexible Endoscopy Order No. Article Description Additional Description Picture

PD-PS-0144p	Modular Catheter 8 F 185cm	for Optic 300 cm, working length: 185 cm integrated light, working channel 1.2mm  disposable / sterile packed packing unit: 5 pcs	
PD-PS-0153p	Modular Catheter 8 F 30cm	for Optic 150 cm, working length: 30 cm integrated light, working channel 1.2mm  disposable / sterile packed packing unit: 5 pcs	0
PD-PS-0154p	Modular Catheter 8 F 42cm	for Optic 150 cm, working length: 42 cm integrated light, working channel 1.2mm  disposable / sterile packed packing unit: 5 pcs	000
PD-PS-0174p	Modular Catheter 8 F 70cm	for Optic 180 cm, working length: 70 cm integrated light, working channel 1.2mm  disposable / sterile packed packing unit: 5 pcs	0-0
PD-TI-6120	tipless basket 70	tipless basket for PolyScope up to 70 cm diameter: 2.5 F length: 120 cm <b>disposable / sterile packed</b>	
PD-TI-6124	4 wire basket 70	4 wire basket for PolyScope up to 70 cm diameter: 2.5 F length: 120 cm disposable / sterile packed	
PD-TI-6214	4 wire basket 185	4 wire basket for PolyScope up to 185 cm diameter: 2.5 F length: 210 cm disposable / sterile packed	
PD-TI-8201	biopsy forceps PolyScope short	biopsy forceps for PolyScope catheter 30 to 70cm diameter: 1.0 mm, length: 120 cm oval spoon  autoclavable	

Please notice: errors excepted; subject and prices to change without notice; similar illustrations





# Flexible Endoscopy Order No. Article Description Additional Description Picture

PD-TI-8202	biopsy forceps PolyScope long	biopsy forceps for PolyScope catheter 185cm diameter: 1.0 mm, length: 220 cm oval spoon autoclavable	
PD-TI-8203	biopsy forceps 0.8 long	biopsy forceps for PolyScope catheter 185cm diameter: 0.8 mm, length: 220 cm oval spoon autoclavable	
PD-TI-9901	oil spray	oil spray for maintenance suitable for all forceps	



Electronics			
Order No.	Article Description	Additional Description	Picture

PD-TR-13953	PolyPort rack 230 Volt	with medically approved power strip for 230 Volt	
PD-VC-0210	camera EndoGnost® CS 207 PAL	high performance CCD camera for all endoscopic applications advanced brightness control with high speed shutter and Moirée filter suitable for 230 Volt / PAL	B 0/1
PD-VC-1022	TV-adaptor 35 mm autoclavable	soakable and autoclavable for easy cleaning	
PD-VR-0401	digital recorder	easy full- featured video recording solution saves videos and images from the endoscope to USB stick (one included) Video in/out: S- video or composite BNC standard 3.5mm mini jack for foot switch (not included) connection please order separately: PD-VR-0402 (foot switch)	2 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0
PD-LS-0220	EndoGnost® LS 200 xenon light source	Xenon light source for a bright, sun- like daylight (5800°K) with brightness fine adjustment xenon lamp with high light yield and life cycle suitable for 230 Volt with IR/UV/Temp. Filter  Suitable for all applications incl. Ophthalmology	Control of the Contro
PD-IP-0100	irrigation pump IP 200	peristaltic pump with adjustable flow rate easy handling and changing of the tubes tubes (PD-IP-0120) and spikes (PD-IP-0121) are available as single use disposables suitable for 230 Volt including foot switch, irrigation bottle with mounting and 5 pcs. PD-IP-0120 and PD-IP-0121 each	



### Catalogue

Electronics			
Order No.	Article Description	Additional Description	Picture

PD-OA-5129	micro drill system complete	micro drill system for semirigid endoscopy recommended espacially for obstruction of lacrimal duct in the tear channel consisting of drive assembly with foot switch, gear transmission piece, driveshaft, drill shifter and motor mount for the 3- joint- arm	
------------	--------------------------------	--	--



# Acessories for Electronics Order No. Article Description Additional Description Picture

PD-VR-0402	footswitch for digital recorder	footswitch for easy handling of the digital recorder	
PD-LC-2001	light cable	length: 180 cm fiber bundle 4.8 mm colour: blue without adaptor	
PD-LC-3001	light cable set	length: 180 cm fiber bundle 4.8 mm colour: blue incl. lightcable PD-LC-2001, adapter Storz PD-LC-9001, adapter ACMI PD-LC-9101	
PD-LC-9001	adaptor C-LS Storz	adaptor for light cable PD-LC-2001 to fit to a Storz light source connection	
PD-LC-9002	adaptor C-LS Fuji	adaptor for light cable PD-LC-2001 to fit to a Fuji light source connection	
PD-LC-9003	adaptor C-LS Pentax	adaptor for light cable PD-LC-2001 to fit to a Pentax light source connection	
PD-LC-9004	adaptor C-LS Olympus/ACMI	adaptor for light cable PD-LC-2001 to fit to a Olympus/ACMI light source connection	
PD-LC-9005	adaptor C-LS Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf light source connection	



# Acessories for Electronics Order No. Article Description Additional Description Picture

PD-LC-9006	adaptor C-LS Stryker	adaptor for light cable PD-LC-2001 to fit to a Stryker light source connection	
PD-LC-9101	adaptor C-E ACMI	adaptor for light cable PD-LC-2001 to fit to a ACMI endoscope connection	=====
PD-LC-9103	adaptor C-E Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf endoscope connection	
PD-LC-9104	adaptor C-E Storz	adaptor for light cable PD-LC-2001 to fit to a Storz endoscope connection	
PD-LC-9502	adaptor E Wolf	adaptor for a PolyDiagnost optic to fit to a Wolf light cable	
PD-LC-9503	adaptor E Storz	adaptor for a PolyDiagnost optic to fit to a Storz light cable	
PD-IP-0120	tube for irrigation pump	tube for irrigation pump PD-IP-0100  disposable / sterile packed	
PD-IP-0121	spike for tube	spike for tube PD-IP-0120  disposable / sterile packed	"₩ <b>()</b>



PD-AC-0101	table clamp	clamps on tables and tubes	THE
PD-AC-0202	vario arm friction	vario arm with knob for variable friction and locking	
PD-AC-0103	camera support	camera support for 3-joint-arm with pins to attach mounting length: 60 cm	-
PD-AC-0118	camera support short	shorter version of PD-AC-0103 length: 50 cm recommended for sales agents	+
PD-AC-0114	mounting	to attach modular oculars and adaptor separate light to the camera support	H
PD-LS-9901	xenon bulb	xenon bulb for light source suitable for PD-LS-0220 and PD-VR-0100	<b>a</b>
PD-IP-0122	foot switch for IP 200	foot switch for irrigation pump PD-IP-0100	00
PD-IP-0123	bottle for irrigation pump	bottle for irrigation pump PD-IP-0100	
PD-OA-5141	drill shifter	holding fixture for the micro drill with Luer Lock connection for adaption to the handle	



### Catalogue

Spare Parts				
Order No.	Article Description	Additional Description	Picture	

PD-OA-5142	driveshaft	flexible spindle for micro drill system PD-OA- 5129 connects gear transmission piece and drill shifter	
PD-OA-5143	gear transmission piece	connects motor of the micro drill system to the driveshaft	
PD-OA-5144	motor mount for 3- Joint-Arm	motor mount for drill system PD-OA-5129 to attach the motor to the 3-Joint-Arm	
PD-VR-1008	flash drive 8GB	PolyDiagnost flash drive for digtal recording suitable for PD-VR-0100 size: 8 GB	