




Simultaneous biplane – efficient and confident prostate examinations


| 8808e |



- Precise biopsies with simultaneous biplane
- Streamlined workflow
- Slender design



Simultaneous biplane



Sagittal plane

Unique patented needle guide

■ Precise biopsies with simultaneous biplane


The 8808e gives you real-time images of both the sagittal and transverse planes, which is invaluable for orientation and precision. Simultaneous images of the sagittal and transverse planes provide a clear indication of correct needle placement for quicker and more confident biopsies.

■ Streamlined workflow

Put efficiency into your routine with the 8808e. Save time and effort with the transducer's one touch operation. Switch between views, freeze, print and save without touching the scanner.

■ Slender design

The 8808e's needle guide is an integral part of the transducer, minimizing patient discomfort.



BK Medical's reusable and sterile, single-use needle guides for 8808e.

Specifications

Frequency range
5 - 10 MHz

Scanning modes
B, M, BCFM, Doppler,
Tissue Harmonic imaging

Image field
126°

Weight
250 g



■ World Headquarters:

Mileparken 34 · DK-2730 · Herlev · Denmark
Tel: (45) 4452 8100 · Fax: (45) 4452 8199 · www.bkmed.com
A wholly owned subsidiary of Analogic Corporation